

ARYS™ pro | starter

Getting better. Everyday.

PRODUCT INFORMATION



What is ARYS™ pro?

Professional therapy support system for the treatment of arm and hand function due to hemiparesis after a stroke or traumatic brain injury. It integrates therapy units into patient's everyday life. The device allows the patient to actively execute and enhance even minimal arm movements which are recorded and contribute to his rehabilitation.

How does ARYS™ pro work?

It is recording arm motion via two wrist-worn trackers. One is recording movements of the healthy arm, the other of the impaired arm. An innovative evaluation software solution running on an iPad tablet transfers the movement records of both trackers into meaningful data, and provides objective, evidence-based documentation for therapists but also for healthcare insurances.

Benefits

The software allows to create single accounts for each patient for individual therapy management, storing patient settings, logging activity, standard tests & results.

Who is ARYS™ pro for?

For healthcare therapists and professionals who treat patients suffering from hemiparesis after a stroke or traumatic brain injury with therapy goals focusing on arm and hand rehabilitation.



We empower rehabilitation.

SYSTEM COMPONENTS



ARYS™ pro | software

available for free on App Store



ARYS™ pro | tracker paretic

2 pcs in ice blue color



ARYS™ pro | tracker healthy

2 pcs in sparkling silver color



ARYS™ pro | band

20 pcs each of them for single patient use



ARYS™ pro | label

20 pcs



ARYS™ pro | tracker charger

Including USB cable & wall hub

ARYS™ pro has been developed in close collaboration with leading neurologists and therapists in the field of neurorehabilitation.

ARYS™ pro has been developed and manufactured according to the applicable standards and the tracker are protected against water splashes.

Manufacturer

yband therapy AG

Laufenstrasse 66

4053 Basel

Switzerland

+41 61 515 64 80

info@ybandtherapy.com